

Test Submission Checklist: Hereditary Breast and Ovarian Cancer BRCA^{now}TM

It's Your Decision

www.GoPathGenetics.com

In order to avoid testing delays, please follow these instructions:

1

Complete the GoPath Hereditary Requisition Form

- Complete family history section
- Is the requisition signed by all required parties?
- Did you obtain copies of insurance cards?
- Do you have all other necessary paperwork?

2

Collect 4ml blood in a purple top EDTA tube



Be sure to write the patient name/ID, date of birth, and specimen collection date on the vial in ink

4

Assemble shipping envelope and ship to GoPath Labs



3

Prepare sample(s) for shipping

Store in refrigerator if specimen is unable to ship the same day of collection. **DO NOT FREEZE**

Inset a refrigerated - NOT FROZEN[®] cold pack into biohazard bag. Place filled tube into bag

Include Hereditary Requisition paperwork, patient's "facesheet"

Place kit in a FedEx Pak and complete shipping label



**Caution:
DO NOT
FREEZE!**

Ship to GoPath Laboratories:
1351 Barclay Boulevard, Buffalo Grove IL 60089

If you do not have a regularly scheduled FedEx pickup, call 800-463-3339 to request one